

Product Datasheet

URM1 antibody (orb166979)

Description	Mouse monoclonal antibody to URM1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	URM1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgM,K
Clonality	Monoclonal
MW	11380
Uniprot ID	Q9BTM9
NCBI	NP_001129419.1 , NP_112176.1
Dilution Range	WB: 1:100
Expiration Date	12 months from date of receipt.

