

Product Datasheet

RYK antibody (orb166993)

Description	Mouse monoclonal antibody to RYK.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	RYK
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2b,k
Clonality	Monoclonal
MW	67815
Uniprot ID	P34925
NCBI	NP_002949.2 , NP_001005861.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

