

Product Datasheet

Human Exosome complex exonuclease RRP44 (DIS3) ELISA Kit (orb1670008)

Description Human Exosome complex exonuclease RRP44 (DIS3) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.