

## Product Datasheet

# Mouse Exosome complex exonuclease RRP44 (DIS3) ELISA Kit (orb1670009)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Mouse Exosome complex exonuclease RRP44 (DIS3) ELISA Kit   |
| <b>Reactivity</b>          | Mouse  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.  |

---

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)