

Product Datasheet

Human Delta sleep-inducing peptide (DSIP) ELISA Kit (orb1670173)

Description	Human Delta sleep-inducing peptide (DSIP) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 18.52-1500 pg/mL. Sensitivity: 7.15 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.