

Product Datasheet

Human Dual Specificity Plaosplaatase (DUSP1) ELISA Kit (orb1670175)

Description	Human Dual Specificity Plaosplaatase (DUSP1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.25-2000 pg/mL. Sensitivity: 12.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.