

Product Datasheet

Human Free Tri-iodothyronine Index (fT3) ELISA Kit (orb1670262)

Description	Human Free Tri-iodothyronine Index (fT3) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 2-32 pmol/L. Sensitivity: 0.8 pmol/L. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.