

Product Datasheet

Rat 4-Hydroxynonenal (HNE) ELISA Kit (orb1670352)

Description	Rat 4-Hydroxynonenal (HNE) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.25-2000 pmol/mL. Sensitivity: 15.5 pmol/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.