

Product Datasheet

Plant Indole-3-acetic acid (IAA) ELISA Kit (orb1670370)

Description	Plant Indole-3-acetic acid (IAA) ELISA Kit
Reactivity	Plant
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.5-800 ng/mL. Sensitivity: 6.5 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.