

Product Datasheet

Mouse Insulin-like growth factor 1 (IGF-1) ELISA Kit (orb1670393)

Description	Mouse Insulin-like growth factor 1 (IGF-1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.25-2000 pg/mL. Sensitivity: 16.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.