

Product Datasheet

ARSA antibody (orb167041)

Description	Mouse monoclonal antibody to ARSA (C-term).
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ARSA
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	53588
Uniprot ID	P15289
NCBI	NP_000478.3 , NP_001078894.2 , NP_001078896.2 , NP_001078895.2 , NP_001078897.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

