

Product Datasheet

Sheep Interleukin-8 (IL-8) ELISA Kit (orb1670465)

Description	Sheep Interleukin-8 (IL-8) ELISA Kit
Reactivity	Sheep
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 7.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.