

Product Datasheet

Human Kallikrein 10 (KLK10) ELISA Kit (orb1670522)

Description	Human Kallikrein 10 (KLK10) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 123.46-10000 pg/mL. Sensitivity: 48.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.