

Product Datasheet

Human Kallikrein-12 (KLK12) ELISA Kit (orb1670524)

Description Human Kallikrein-12 (KLK12) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 12.35-1000 ng/mL.
Sensitivity: 7.14 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.