



Product Datasheet

Plp1 antibody (orb167056)

Description	Mouse monoclonal antibody to Plp1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Plp1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	30077
Uniprot ID	P60203
NCBI	NP_112252.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse cerebellu...