

Product Datasheet

Human Extracellular glycoprotein lacritin (LACRT) ELISA Kit (orb1670562)

Description	Human Extracellular glycoprotein lacritin (LACRT) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.63 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.