

Product Datasheet

Mouse L-lactate dehydrogenase B chain (LDHB) ELISA Kit (orb1670592)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Mouse L-lactate dehydrogenase B chain (LDHB) ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.156-10 ng/mL. Sensitivity: 0.055 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids. |
| Expiration Date | Please enquire. |