

Product Datasheet

Swine IL-2 ELISA KIT (orb1670746)

Description	Swine IL-2 ELISA KIT
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.02-1.5 ng/mL. Sensitivity: 12 pg/mL. Assay Time: 3 hour. Sample Volume: 20 μ L. Sample Type: Cell culture supernatant, serum, plasma (EDTA, citrate, heparin). Detection method: Colorimetric.
Expiration Date	Please enquire.