

Product Datasheet

Human PBEF/ Visfatin ELISA KIT (orb1670784)

Description	Human PBEF/ Visfatin ELISA KIT
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 pg/mL. Sensitivity: 0.70 ng/mL. Assay Time: 3 hours. Sample Volume: 20 µL. Sample Type: Cell culture supernatant, serum, plasma (EDTA, citrate, heparin). Detection method: Colorimetric.
Expiration Date	Please enquire.