

Product Datasheet

Sheep Follicle Stimulating Hormone (FSH) ELISA kit (orb1670827)

Description	Sheep Follicle Stimulating Hormone (FSH) ELISA kit
Reactivity	Sheep
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0 mIU/mL-70 mIU/mL. Sensitivity: 1 mIU/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.