

## Product Datasheet

### Rat Tri-iodothyronine, T3 ELISA Kit (orb1670837)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rat Tri-iodothyronine, T3 ELISA Kit  |
| <b>Reactivity</b>          | Rat  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detect Range: 0.5 ng/mL-8 ng/mL.<br>Sensitivity: 0.25 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm.<br>Sample Volume: 50-100 $\mu$ L. Sample Type: serum, plasma, cell culture supernates, tissue homogenates. |
| <b>Expiration Date</b>     | Please enquire.  |