

Product Datasheet

Rat Progesterone,PROG ELISA kit (orb1670839)

Description	Rat Progesterone,PROG ELISA kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.2 ng/mL-50 ng/mL. Sensitivity: 0.2 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, cell culture supernates, tissue homogenates.
Expiration Date	Please enquire.