



Product Datasheet


Myc Tag antibody (orb167084)

Description	Mouse monoclonal antibody to Myc Tag.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	KLH conjugated synthetic peptide encoding c-Myc tag (AEEQKLISEEDLLRKRREQLKHKLE) was used as antigen.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell line

...



Western blot analysis of mouse Neuro-2a ...