

Product Datasheet

Rat microalbuminuria(MAU/ALB) ELISA kit (orb1670840)

Description	Rat microalbuminuria(MAU/ALB) ELISA kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.01 µg/mL-40 µg/mL. Sensitivity: 0.01 µg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, urine.
Expiration Date	Please enquire.