

Product Datasheet

Rat Free Tri-iodothyronine Index, Free-T3 ELISA Kit (orb1670845)

Description	Rat Free Tri-iodothyronine Index, Free-T3 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 2 pmol/L-32 pmol/L. Sensitivity: 0.38 pmol/L. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates.
Expiration Date	Please enquire.