

Product Datasheet

Rat Free Thyroxine, FT4 ELISA Kit (orb1670846)

Description	Rat Free Thyroxine, FT4 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 4 pmol/L-64 pmol/L. Sensitivity: 1.08 pmol/L. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma.
Expiration Date	Please enquire.