



Product Datasheet

HA Tag antibody (orb167085)

Description	Mouse monoclonal antibody to HA Tag.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	KLH conjugated synthetic peptide encoding HA tag (YPYDVPDYA) was used as antigen.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 12tag recombina...