

Product Datasheet

Rabbit Albumin(Alb) ELISA Kit (orb1670866)

Description	Rabbit Albumin(Alb) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.273 µg/mL-70 µg/mL. Sensitivity: 0.117 µg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, urine.
Expiration Date	Please enquire.