

Product Datasheet

Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit (orb1670869)

Description	Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.047 ng/mL-3 ng/mL. Sensitivity: 0.040 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, cell culture supernates, tissue homogenates.
Expiration Date	Please enquire.