




Product Datasheet

HIS Tag antibody (orb167087)

Description	Mouse monoclonal antibody to HIS Tag.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	Purified recombinant HIS-tagged fusion protein and poly-HIS peptide were used to produced this monoclonal antibody.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis of 12 Tag protein using HIS ...



WB analysis of 12 Tag protein using HIS ...