



Product Datasheet

GST antibody (orb167089)

Description	Mouse monoclonal antibody to GST.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	Purified recombinant GST fusion protein was used to produced this monoclonal antibody.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1 k
Clonality	Monoclonal
MW	26000
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of GST recombinant...