

Product Datasheet

Mouse Pentraxin-related protein PTX3(PTX3) ELISA kit (orb1670897)

Description	Mouse Pentraxin-related protein PTX3(PTX3) ELISA kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.313 ng/mL-20 ng/mL. Sensitivity: 0.210 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.