

Product Datasheet

Mouse hepatitis B virus e antibody, HBeAb ELISA Kit (orb1670907)

| | |
|----------------------------|---|
| Description | Mouse hepatitis B virus e antibody, HBeAb ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma. |
| Expiration Date | Please enquire. |