

## Product Datasheet

# Mouse hepatitis B virus core antibody, HBcAb ELISA Kit (orb1670908)

<b>Description</b>	Mouse hepatitis B virus core antibody, HBcAb ELISA Kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 $\mu$ L. Sample Type: serum, plasma.
<b>Expiration Date</b>	Please enquire.