

## Product Datasheet

# Mouse anti respiratory syncytial virus (RSV) antibody (IgG)ELISA kit (orb1670924)

<b>Description</b>	Mouse anti respiratory syncytial virus (RSV) antibody (IgG)ELISA kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 $\mu$ L. Sample Type: serum.
<b>Expiration Date</b>	Please enquire.