

Product Datasheet

Human Vascular endothelial cell growth factor receptor 2, VEGFR-2/FIk-1 ELISA kit (orb1670941)

Description	Human Vascular endothelial cell growth factor receptor 2, VEGFR-2/FIk-1 ELISA kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.39 ng/mL-25 ng/mL. Sensitivity: 0.098 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.