

Product Datasheet

Human soluble cluster of differentiation 14,sCD14 ELISA Kit (orb1670958)

Description Human soluble cluster of differentiation 14,sCD14 ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.234 ng/mL-15 ng/mL. Sensitivity: 0.229 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, cell culture supernates, tissue homogenates and urine.

Expiration Date Please enquire.