

Product Datasheet

Human microalbuminuria(MAU/ALB) ELISA kit (orb1670976)

Description	Human microalbuminuria(MAU/ALB) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1 µg/mL-50 µg/mL. Sensitivity: 0.5 µg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, urine.
Expiration Date	Please enquire.