

Product Datasheet

Human Macrophage-Derived Chemokine (MDC/CCL22) ELISA kit (orb1670982)

Description	Human Macrophage-Derived Chemokine (MDC/CCL22) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.625 pg/mL-1000 pg/mL. Sensitivity: 10.324 pg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, cell culture supernates, tissue homogenates.
Expiration Date	Please enquire.