

Product Datasheet

Human Folic acid(FA) ELISA Kit (orb1671013)

Description	Human Folic acid(FA) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.5 pg/mL-200 pg/mL. Sensitivity: 6.25 pg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma.
Expiration Date	Please enquire.