

Product Datasheet

Human Epithelial neutrophil activating peptide 78 (ENA-78/CXCL5) ELISA Kit (orb1671017)

Description	Human Epithelial neutrophil activating peptide 78 (ENA-78/CXCL5) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625 ng/mL-40 ng/mL. Sensitivity: 0.567 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.