

Product Datasheet

Human Cystatin-B (CSTB/CST6/STFB) ELISA kit (orb1671023)

Description	Human Cystatin-B (CSTB/CST6/STFB) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.313 ng/mL-20 ng/mL. Sensitivity: 0.233 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, cell culture supernates, urine, saliva, tissue homogenates, cell lysates.
Expiration Date	Please enquire.