

Product Datasheet

Human Chitinase-3-like Protein 1(YKL-40/CHI3L1)ELISA Kit (orb1671033)

Description	Human Chitinase-3-like Protein 1(YKL-40/CHI3L1)ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 46.875 pg/mL-3000 pg/mL. Sensitivity: 78.065 pg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.