

## Product Datasheet

### Goat Albumin (Alb)ELISA Kit (orb1671114)

<b>Description</b>	Goat Albumin (Alb)ELISA Kit
<b>Reactivity</b>	Goat
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.17 µg/mL-350 µg/mL. Sensitivity: 0.17 µg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, tissue homogenates.
<b>Expiration Date</b>	Please enquire.