

Product Datasheet

Fish free thyroxine (FT4) ELISA kit (orb1671119)

Description	Fish free thyroxine (FT4) ELISA kit
Reactivity	Fish
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 4 pmol/l-64 pmol/l. Sensitivity: 1.08 pmol/L. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma.
Expiration Date	Please enquire.