

Product Datasheet

Canine thyroxine, T4 ELISA Kit (orb1671147)

Description	Canine thyroxine, T4 ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 20 ng/mL-320 ng/mL. Sensitivity: 10 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, cell culture supernates, tissue homogenates.
Expiration Date	Please enquire.