


Product Datasheet

Interleukin 4 Recombinant Protein (orb1671201)

Description	Interleukin 4 Recombinant Protein
Tested Applications	CLIA, ELISA, LF, WB
Preservatives	PBS, pH 8.0, and 8 M Urea containing 0.1% ProClin300.
Form/Appearance	liquid
Concentration	batch dependent
Storage	Upon receipt, store at 2-8°C
Note	For research use only
Application notes	Suitable for use in ELISA, CLIA, Lateral Flow, Western Blot, and other immune tests
Purity	> 95% pure by SDS-PAGE. Ligand affinity purified.
Source	E. coli
MW	This recombinant protein contains full length of Human interleukin 4 with a MW of 16.2 kDa. It is recognized by an anti-human interleukin 4 mAb in ELISA tests.
Expiration Date	6 months from date of receipt.



SDS-Page
Image