


## Product Datasheet

# Interferon (IFN) gamma Recombinant Protein (orb1671207)

|                            |                                                                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Interferon (IFN) gamma Recombinant Protein                                                                                                            |
| <b>Tested Applications</b> | CLIA, ELISA, LF, WB                                                                                                                                   |
| <b>Preservatives</b>       | PBS, pH 8.0, and 8 M Urea containing 0.1% ProClin300.                                                                                                 |
| <b>Form/Appearance</b>     | liquid                                                                                                                                                |
| <b>Concentration</b>       | batch dependent                                                                                                                                       |
| <b>Storage</b>             | Upon receipt, store at 2-8°C                                                                                                                          |
| <b>Note</b>                | For research use only                                                                                                                                 |
| <b>Application notes</b>   | Suitable for use in ELISA, CLIA, Lateral Flow, Western Blot, and other immune tests                                                                   |
| <b>Purity</b>              | > 90% pure by SDS-PAGE. Ligand affinity purified.                                                                                                     |
| <b>Source</b>              | E. coli                                                                                                                                               |
| <b>MW</b>                  | This recombinant protein contains full length of Human IFN-gamma with a MW of 18 kDa. It is recognized by an anti-human IFN-gamma mAb in ELISA tests. |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                        |



SDS-Page  
Image