


Product Datasheet

Anti-Müllerian Hormone (AMH) human Recombinant Protein (orb1671218)

Description	Anti-Müllerian Hormone (AMH) human Recombinant Protein
Tested Applications	CLIA, ELISA, LF, WB
Preservatives	PBS, pH 7.3 containing 0.1% ProClin300.
Form/Appearance	liquid
Concentration	batch dependent
Storage	Upon receipt, quick spin before open vial, aliquot and store at $\leq -20^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles
Note	For research use only
Application notes	Suitable for use in ELISA, CLIA, Lateral Flow, Western Blot, and other immune tests
Purity	> 90% pure by SDS-PAGE. Ligand affinity purified.
Source	Insect cell
MW	This recombinant protein contains full length of human AMH with a calculated MW of 61 kDa. It is recognized by anti-human AMH antibodies from R&D systems in ELISA tests.
Expiration Date	6 months from date of receipt.



SDS-
Page
Image