

## Product Datasheet

# Plasmodium falciparum histidine-rich 2 (Pf-HRP II) Recombinant Protein (orb1671332)

**Description** Plasmodium falciparum histidine-rich 2 (Pf-HRP II) Recombinant Protein

**Tested Applications** CLIA, ELISA, LF, WB

**Preservatives** Tris-buffered saline, pH 8.0 containing 0.1% ProClin300.

**Form/Appearance** liquid

**Concentration** batch dependent

**Storage** Upon receipt, quick spin before open vial, aliquot and store at  $\leq -20^{\circ}\text{C}$ . Avoid multiple freeze/thaw cycles

**Note** For research use only


**Application notes** Suitable for use in ELISA, CLIA, Lateral Flow, Western Blot, and other immune tests

**Purity** > 90% pure by SDS-PAGE. Ligand affinity purified.

**Source** E. coli

**MW** This recombinant protein contains full length of Pf-HRP II with a MW of 56.4 kDa. It is recognized by positive samples in ELISA tests.

**Expiration Date** 6 months from date of receipt.



SDS-  
Page  
Image